Liquid Anti-Xa - 0020302600

Intended use

Automated chromogenic assay for the quantitative determination of unfractionated heparin (UFH) and low molecular weight heparin (LMWH) activity in human citrated plasma on IL Coagulation Systems, (ACL TOP® Family, ACL TOP Family 50 Series¹ and ACL Elite®/Elite Pro/8/9/10000).

Summary and principle

Heparin is the most frequently used antithrombotic drug. The biological activity of this sulphated glycosaminoglycan resides in its ability to accelerate (up to 2000-fold) the inhibitory effect of antithrombin on coagulation proteases. In recent years, it has been shown that LMWH, besides being as useful therapeutically as UFH, also has a longer half-life.

The Liquid Anti-Xa kit is a one stage chromogenic assay based on a synthetic chromogenic substrate and on Factor Xa inactivation. Heparin levels in patient plasma are measured automatically on IL Coagulation Systems.

Heparin is analyzed as a complex with antithrombin present in the sample. The concentration of this complex is dependent on the availability of the patient's endogenous antithrombin. When the Heparin – antithrombin complex is formed, two competing reactions take place.

- Factor Xa is neutralized by heparin-antithrombin complex.
- Residual Factor Xa is quantified with a synthetic chromogenic substrate. The paranitroaniline released is monitored kinetically at 405 nm and is inversely proportional to the heparin level in the sample.¹

In order to reduce the influence from heparin antagonists, such as platelet factor 4 (PF4), dextran sulfate is included in the reaction mixture.

Composition

The **Anti-Xa** kit consists of:

[E] **Factor Xa reagent** (Cat. No. 0020302610): 5 x 2.5 mL vial of a liquid preparation containing purified bovine Factor Xa (approximately 5.5 nkat/mL), Tris-Buffer, EDTA, dextran sulfate, sodium chloride and bovine serum albumin.

[S] **Chromogenic substrate** (Cat. No. 0020302620): 5 x 3 mL vial of liquid chromogenic substrate S-2732 (approximately 1.2 mg/mL) and bulking agent.

PRECAUTIONS AND WARNINGS:

The FXa Reagent contains bovine material. All donor animals were sourced from BSE-free herds. The cattle received ante- and post mortem health inspection by a veterinarian, and they were apparently free from infectious and contagious material. However, the material should be treated as potentially infectious.

Factor Xa reagent

Danger
Hazard class: Resp. Sens. 1, H334
Hazard statements: H334: May cause allergy or asthma symptoms or breathing difficulties if inhaled.

Precautionary statements: P261: Avoid breathing vapors/ spray. P284: [In case of inadequate ventilation] wear respiratory protection. P304 + P340: IF INHALED: Remove person to fresh air and keep comfortable for breathing. P342 + P311: If experiencing respiratory symptoms: Call a POISON CENTER/doctor. P501: Dispose of contents/container in accordance with local/regional/national/international regulation.

Supplemental hazard information: Contains Factor Xa. Up to 2.8% of the mixture consists of component of unknown acute toxicity (oral, dermal, inhalation) for the human health and for the aquatic environment.**Chromogenic substrate**

Hazard class: none

Hazard statements: none

Precautionary statements: none

Supplemental hazard information: none

This product is **For in vitro** Diagnostic Use.

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Preparation

Chromogenic substrate: Invert to mix before use.

Factor Xa reagent: Invert to mix before use.

Reagent storage and stability

Unopened reagents are stable until the expiration date shown on the vial when stored at 2-8°C.

Chromogenic substrate – Opened reagent is stable: 1 month at 2-8°C or 7 days at 15° -25°C on-board the ACL TOP Family and ACL TOP Family 50 Series¹, 3 days on-board the ACL Elite/Elite Pro/8/9/10000 Systems in the original vial.

Factor Xa reagent - Opened reagent is stable: 1 month at 2-8°C or 7 days at 15° -25°C on-board the ACL TOP Family and ACL TOP Family 50 Series¹, 3 days on-board the ACL Elite/Elite Pro/8/9/10000 Systems in the original vial.

For optimal stability remove reagents from the system and store them at 2-8°C in the original vial.

⚠ Low relative humidity is associated with increased evaporation of uncapped reagents, which may decrease on-board stability. For optimum on-board stability, laboratory temperature and humidity should be controlled.

Instrument/test procedures

Refer to the appropriate IL instrument's Operator's Manual for the complete assay procedure instructions.

Specimen collection and preparation

Nine parts of freshly drawn venous blood are collected into one part 3.2% trisodium citrate. Refer to CLSI Document H21-A5 for further instructions on specimen collection, handling and storage.²

Additional reagent and control plasmas

The following are not supplied with the kit and must be purchased separately.

	Product Cat No.
Heparin Calibrators	0020300600
LMW Heparin Controls	0020300200
UF Heparin Controls	0020300300
Cleaning solution	0009831700
Cleaning agent	0009832700

Quality control

Two levels of controls are recommended for a complete quality control program.³ LMW and UF Heparin Controls are designed for this program. Each laboratory should establish its own mean and standard deviation and should establish a quality control program to monitor laboratory testing. Controls should be analyzed at least once every 8 hour shift in accordance with good laboratory practice. Refer to the instrument's Operator's Manual for additional information. Refer to Westgard *et al*/ for identification and resolution for out-of-control situations.⁴

Results

Heparin results are reported in IU/mL. Refer to the instrument's Operator's Manual for additional information.

Limitations/interfering substances

Heparin results on the ACL TOP Family, ACL TOP Family 50 Series¹ and ACL Elite/Elite Pro/8/9/10000 Systems are not affected by hemoglobin up to 300 mg/dL, bilirubin up to 20 mg/dL and triglycerides up to 800 mg/dL.

Expected values

To obtain an optimal effect with minimum risk of bleeding or thromboembolic complications the heparin activity should be in the range recommended by the heparin manufacturer.⁵

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ENGLISH - Insert revision 06/2017

Performance characteristics

Precision
Within run and total (run - to - run and day - to - day) precision was assessed over multiple runs using the two levels of both UFH and LMWH Controls.

ACL TOP Family	Mean (IU/mL)	CV% (Within run)	CV% (Total)
UFH Low	0.41	2.4	3.3
UFH High	0.68	1.0	1.6
LMWH Low	0.55	3.5	4.5
LMWH High	1.35	1.9	2.5
ACL TOP Family 50 Series¹	Mean (IU/mL)	CV% (Within run)	CV% (Total)
UFH Low	0.33	1.9	4.0
UFH High	0.63	1.3	2.0
LMWH Low	0.62	1.5	3.6
LMWH High	1.51	1.2	2.5
ACL Elite/Elite Pro/8/9/10000	Mean (IU/mL)	CV% (Within run)	CV% (Total)
UFH Low	0.39	4.0	6.7
UFH High	0.69	1.0	3.7
LMWH Low	0.49	5.5	5.7
LMWH High	1.31	3.3	4.0

Correlation: System	Slope	Intercept	r	Reference method
ACL TOP Family	0.946	0.055	0.958	HemosIL Heparin
ACL TOP Family 50 Series¹	0.973	-0.020	0.999	HemosIL Liquid Anti-Xa on ACL TOP 500 CTS HemosIL Heparin
ACL Elite/Elite Pro/8/9/10000	0.894	0.057	0.907	

The precision and correlation results were obtained using specific lots of reagent and controls. In field site studies using samples from patients undergoing heparin therapy, the following data were obtained using a specific lot of Liquid Anti-Xa reagents:

System	Slope	Intercept	r	Reference method
ACL TOP Family	0.952	0.095	0.978	HemosIL Heparin
ACL Elite/Elite Pro/8/9/10000	1.032	0.039	0.949	HemosIL Heparin

Detection limit: System

ACL TOP Family/ACL TOP Family 50 Series¹ 0.04 IU/mL
ACL Elite/Elite Pro/8/9/10000 0.04 IU/mL

Linearity: System

ACL TOP Family/ACL TOP Family 50 Series¹ Up to - 2.0 IU/mL
ACL Elite/Elite Pro/8/9/10000 Up to - 2.0 IU/mL

† ACL TOP Family 50 Series = ACL TOP 350 CTS; ACL TOP 550 CTS; ACL TOP 750; ACL TOP 750 CTS; ACL TOP 750 LAS

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- Tieen AN, *et al*. Assay of Heparin in Plasma using a Chromogenic Substrate, Thromb. Res. 1976; 8: 413-416.
- Clinical and Laboratory Standards Institute/CLSI. Collection, Transport and Processing of Blood Specimens for Testing Plasma-Based Coagulation and Molecular Haemostasis Assays; Approved Guideline – Fifth Edition, CLSI Document H21-A5; Vol. 28 No. 5.
- Zucker S, Cathey MH, West B. Preparation of Quality Control Specimens for Coagulation, Am J. Clin. Pathol. 1970; 53: 924-927.

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