Bestellinformation

REF	CONTENT	Gerät(e), auf dem/denen das cobas c pack/die cobas c packs verwendet werden kann/können
06600239 190	Tina-quant Cystatin C Gen.2 (225 Tests)	System-ID 07 7550 9 Roche/Hitachi cobas c 311, cobas c 501/502
04975901 191	C.f.a.s. Cystatin C (4 x 1 mL)	Code 407
04975936 190*	Cystatin C Control Set Kontrolle I (niedrig) (4 x 1 mL) Kontrolle II (hoch) (4 x 1 mL)	Code 121 Code 122
06729371 190	Cystatin C Control Set Gen.2 Kontrolle 1 (3 x 1 mL) Kontrolle 2 (3 x 1 mL) Kontrolle 3 (3 x 1 mL)	Code 139 Code 140 Code 141

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Deutsch

Systeminformation

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Literatur

- Levey AS, Coresh J, Balk E, et al. National Kidney Foundation practice guidelines for chronic kidney disease: evaluation, classification, and stratification. *Ann Intern Med* 2003;139:137-147.
- Rule AD, Larson TS, Bergstralh EJ, et al. Using serum creatinine to estimate glomerular filtration rate: accuracy in good health and in chronic kidney disease. *Ann Intern Med* 2004;141:929-937.
- Wasen E, Isoaho R, Mattila K, et al. Estimation of glomerular filtration rate in the elderly: a comparison of creatinine based formulae with serum Cystatin C. *J Intern Med* 2004;256:70-78.
- Levey AS, Coresh JK. DOQI clinical practice guidelines on chronic kidney disease. Guideline 4. Estimation of GFR. *Am J Kidney Dis* 2002;39 (suppl1):76-92.
- Kyhse-Andersen J, Schmidt C, Nordin G, et al. Serum Cystatin C, determined by a rapid, automated particle-enhanced turbidimetric method, is a better marker than serum creatinine for glomerular filtration rate. *Clin Chem* 1994;40(10):1921-1926.
- Grubb A, Simonsen O, Sturfelt G, et al. Serum concentration of Cystatin C, factor D and β 2-microglobulin as a measure of glomerular filtration rate. *Acta Med Scand* 1985;218(5):499-503.
- Jung K, Jung M. Cystatin C: a promising marker of glomerular filtration rate to replace creatinine. *Nephron* 1995;70(3):370-371.
- Newman DJ, Thakkar H, Edwards RG, et al. Serum Cystatin C measured by automated immunoassay: a more sensitive marker of changes in GFR than serum creatinine. *Kidney Int* 1995;47(1):312-318.
- Mussap M, Dalla Vestra M, Fioretto P, et al. Cystatin C is a more sensitive marker than creatinine for the estimation of GFR in type 2 diabetic patients. *Kidney Int.* 2002;61:1453-1461.
- Dhanidharka VR, Kwon C, Stevens G. Serum cystatin C is superior to serum creatinine as a marker of kidney function: a meta analysis. *Am J Kidney Dis* 2002 Aug;40(2):221-226.
- Risch L, Drexel H, Huber AR. Differences in glomerular filtration rate estimates by two cystatin C-based equations. *Clinical Chemistry* 2005;51:2211-2212.
- Grubb A, Nyman U, Bjork J, et al. Simple Cystatin C-Based prediction equations for glomerular filtration rate compared with the modification of diet in renal disease prediction equation for adults and the Schwartz and the Counahan-Barratt prediction equations for children. *Clinical Chemistry* 2005;51:1420-1431.
- Hoek FJ, Kemperman FA, Krediet RT. A comparison between Cystatin C, plasma creatinine and the Cockcroft and Gault formula for the estimation of glomerular filtration rate. *Nephrol Dial Transplant* 2003;18:2034-2031.
- Lassus J, Harjola VP, Sund R, et al. Prognostic value of Cystatin C in acute heart failure in relation to other markers of renal function and NT-proBNP. *European Heart Journal* (February) doi:10.1093/eurheart/ehl507.
- Stevens H, Wide-Svensson D, Torffvit O, et al. Serum Cystatin C for assessment of glomerular filtration rate in pregnant and non-pregnant women. Indications of altered filtration process in pregnancy. *Scand J Clin Lab Invest* 2002;62:141-148.
- Andersen TB, Eskild-Jensen A, Frokiaer J, et al. Measuring glomerular filtration rate in children; can Cystatin C replace established methods? *Pediatr Nephrol* 2009;24:929-41.
- Grubb A, Nyman U, Björk, J Improved estimation of glomerular filtration rate (GFR) by comparison of eGFR (cystatin C) and eGFR (creatinine). *Scand J Clin Lab Invest* 2012;72:73-77.
- Filler G, Huang S-HS, Yasin A.. The usefulness of Cystatin C and related formulae in pediatrics. *Clin Chem Lab Med* 2012;50(12):2081-2091.
- Grubb A. Cystatin C- and creatinine-based GFR-prediction equations for children and adults. *Clinical Biochemistry* 2011;44:501-502.
- Rollins G. Predicting risk of ESRD in Diabetics, *Clinical Laboratory Strategies*, Jan.2013,1-4.
- Leseley A, Inker M.D., Christopher H, et al. Estimating Glomerular Filtration Rate from Serum Creatinine and Cystatin C. *N Engl J Med* 2012;367:20-29.
- Horio M, Imai E, Yasuda Y, et al. GFR Estimation Using Standardized Serum Cystatin C in Japan. *Am J Kidney Dis* 2013;61(2):197-203
- Grubb A, Horio M, Hansson LO, et al. Generation of a New Cystatin C-Based Estimating Equation for Glomerular Filtration Rate Using Seven Assays Standardized to the International Calibrator. *Clin Chem* 2014, in press.
- van Deutekom AW, Zur B, van Wijk JAE, et al. Measurement of cystatin C in capillary blood samples in pediatric patients. *Clin Biochem* 2010;43:335-337.
- Data on file at Roche Diagnostics.
- Gislefoss RE, Grimsrud TK, Morkrid L. Stability of selected serum proteins after long-term storage in the Janus Serum Bank. *Clin Chem Lab Med* 2009;47(5):596-603.
- Bökenkamp A, van Wijk JAE, Lentze ML, et al. Effect of Corticosteroid Therapy on Serum Cystatin C and β 2-Microglobulin Concentrations. *Clin Chem* 2002;48:1123-1126.
- Wiesli P, Schwegler B, Spinass AS, et al. Serum Cystatin C is sensitive to small changes in thyroid function *Clin Chim Acta* 2003;338(1):87-90.
- Glick MR, Ryder KW, Jackson SA. Graphical Comparisons of Interferences in Clinical Chemistry Instrumentation. *Clin Chem* 1986;32:470-475.
- Breuer J. Report on the Symposium "Drug effects in Clinical Chemistry Methods". *Eur J Clin Chem Clin Biochem* 1996;34:385-386.
- Sonntag O, Scholer A. Drug interference in clinical chemistry: recommendation of drugs and their concentrations to be used in drug interference studies. *Ann Clin Biochem* 2001;38:376-385.
- Bakker AJ, Mücke M. Gammopathy interference in clinical chemistry assays: mechanisms, detection and prevention. *Clin Chem Lab Med* 2007;45(9):1240-1243.
- Kricka LJ. Human Anti-Animal Antibody Interferences in Immunological Assays. *Clin Chem* 1999;45(7):942-956.
- Bablok W, Passing H, Bender R, et al. A general regression procedure for method transformation. Application of linear regression procedures for method comparison studies in clinical chemistry, Part III. *J Clin Chem Clin Biochem* 1988 Nov;26(11):783-790.

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© 2016, Roche Diagnostics
Roche Diagnostics GmbH, Sandhofer Strasse 116, D-68305 Mannheim